DIGITAL MICROSCOPE

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Beaverlab

Please read this instruction manual in detail.
Please keep this instruction manual in a safe place.

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
01/ Product Overview

Thank you for purchasing the Digital Microscope, hereinafter referred to as the microscope. Please read this manual carefully before use and keep it in a safe place. Please do not use the microscope without understanding how to use it.


Instructions For Use


- Please read this manual carefully and use it correctly.
- Pay attention to the protection of the lens. Be sure to close the lens cap after use to avoid dust or foreign objects from getting inside.
- This microscope can be used handheld or inserted into the base for desktop use; pay attention to the strength when inserted to avoid damage to the components.
- Place the specimen directly below the center of the lens to avoid not being able to see the object and affecting the observation experience.
- Take care to protect the screen when using. Do not scratch or deface the screen. Clean the screen after each use.
- Import of external files is not supported (external files such as videos, pictures, and documents do not support import and reading, which may cause the device to malfunction).

02/The Microscope Instructions

How To Open And Use The Microscope

Before using the microscope, please make sure to hold it firmly in your hand or place it on a tabletop for use.

Microscope Power Button

Press and hold for 3 seconds to turn on the phone.
Press and hold for 3 seconds to turn off the phone.

Digital Zoom Button

Press briefly to digitally enlarge.

Digital Reduction Key

Press briefly to digitally reduce.

Microscope Monomer Interface

Charge and connect to the computer.

Brightness Dimming Key

Press briefly to brighten.

Brightness Dimmer Key

Press briefly to dim.

Microscope On And Off


- Press and hold the microscope rear power on button for 3 seconds. When the lower small light of the button lights up or the bottom illumination lights up, that means the microscope has been turned on.
- Press and hold the on button at the back of the microscope for 3 seconds. The small light at the bottom of the button will go out or the bottom illumination will go out, which means the microscope has been turned off.

(Note: After using the microscope, please put the lens cap on.)

Microscope Magnification Adjustment

-  Counterclockwise rotation: Zoom in
-  Clockwise rotation: Zoom out

(Note: At the same distance, there are two magnifications by rotating the focusing rotary drum. The lifting knob and digital zoom-in and zoom-out keys can also be used to adjust the magnification of the microscope for observation. The focusing rotary drum, lifting knob, and digital zoom-in function are used together.)

Usage Environment

Please use in -10°C~45°C environment.

Please avoid getting the Microscope wet with water, drinking water, corrosive liquid, etc.

Please avoid placing the microscope near heat, open flames, and flammable and explosive gases/liquids.


Do not let dust enter the microscope lens and interior.

Do not hit the microscope. Avoid violently vibrate the Microscope.

- At low temperatures, the available battery capacity decays to varying degrees.

Charging Instructions

Indicator light

The microscope can be charged with a cell phone charger, computer, or rechargeable battery.

Do not charge for more than 12 hours to avoid affecting the battery life.

- When the microscope battery is low, the power indicator's red light flashes fast.
- The red light is always on when charging, and the light goes off when full.
- During the charging process, it is normal for the microscope case and the battery to heat up slightly, so please don't worry about the warmth.

Storage Environment


- Store in a cool, dry place for daily use and avoid direct sunlight.
- Avoid storing the microscope in places where there is a risk of dropping, which may cause damage to the internal lens deviation or other components of the microscope, as well as other irreparable damage.

Step-By-Step Instructions For Use


1. Hold the lower part of the microscope with one hand.


2. Press and hold the rear power-on button until the screen lights up.


3. Remove the lens cover.

4. Aim the microscope at the object to be observed.

5. The other hand rotates the focus rotating cylinder. Observe the screen display, and adjust to the best observation effect.

03/Accessories and tools

Data Cable


Instruction Manual

04/Interface Usage Instructions



- **Turn On The Device**

Press and hold the power on button for 3 seconds to turn on the phone


- **Photo And Video**

Touch the photo icon to take a picture, and touch the record icon to start/end recording.

• Album

Press the menu button to bring up the three buttons of video album/photo album/settings. Select left and right and press the confirmation button to enter the corresponding interface.

Video Album



Photo Album

Settings

- Date/Time
- Language
- Sleep (Off, 3 min, 5 min, 10 min, 15 min) *Press the button during sleep to return to normal
- Auto-off timer (Off, 5 min, 10 min, 30 min, 60 min)
- Format
- Initialize
- Power off
- WiFi switch (Off/On) *When turned on, you can connect to your smartphone (Not all smartphones can be connected)

Settings (microscope
-related settings)

App Interface User Guide

• App download and installation:

Scan the QR code below to download the app (search "Beaver Point" in App Store or Google Play to download the app) Install the software according to the operating instructions.

Android

iOS

(Scan here to download the app)

• Turn On The Device

Press and hold the microscope rear power on button for 3 seconds. The lower small light of the button lights up or the bottom illumination lights up, which means the microscope has been turned on.

• Connect microscope to Android/iOS device (WIFI initial password: 12345678)

After the device is turned on, open the mobile app. Click to connect to WIFI and find the hotspot of "DDL-M2-XXXXXX." Tap and connect to the hotspot. If a connection confirmation prompt pops up, please click "Stay connected" (Android version) /click "OK." After the connection is successful, the return software will display that the device has been successfully connected, which proves that the microscope has successfully connected with your device.

Click to
connect to WIFI


Select the corresponding
WIFI name


Device successfully
connected

- After clicking Live View, you can enter the real-time observation interface.
(Note: After the app is closed, you need to manually disconnect the WIFI to continue the observation.)

Click to watch live

Real-time observation


For Computer Use

Log on to www.beaverlabtech.com to download and install the software installation package. Install software according to operating instructions.

Use the USB cable that comes with the microscope to connect the microscope tube to the computer host.

Open the software to select microscope equipment.

Follow the software operating instructions to use it. (Tips: After the device is connected to the computer, the software automatically recognises the device and selects the resolution.)

Computer users can go to the official website www.beaverlabtech.com to download the Windows version of the software

05/Troubleshooting Methods

Problem	Reason	Processing Method
Heat during charging	High current during charging, high battery heat, which is normal	Normal state
The app can't connect	①Distance too far ②Equipment without power ③Equipment failure	①Repeat connection in close proximity ②Use after charging the device ③App and device are closed and restarted before use ④Need maintenance
App suddenly disconnected	①Distance too far ②Equipment without power ③Equipment failure	①Repeat connection in close proximity ②Use after charging the device ③App and device are closed and restarted before use ④Need maintenance
Device cannot be charged	①Charging port is not plugged firmly ②No power access ③Equipment failure	①Replug the interface ②Connect the power ③Need maintenance
The device is not working properly	Internal damage caused by water, collision, etc.	Need maintenance
The device cannot connect to the software side of the computer	①The interface is not inserted tightly ②The data cable is damaged ③Not using the connection cable with data transfer function	①Replug the interface ②Replace the data cable ③Replace the data cable with data transfer function and then connect to use
Lens display image blur	Lens fogging due to radical changes in ambient temperature	Wait for the temperature to equilibrate and the fog to disappear before use
Drop, impact resulting in broken and cracked shell	Damage caused by external use	Need maintenance
The product cannot be turned off or the button does not work	Crash	Perform a reset operation
Black screen does not display, flashing screen, etc. .	Malfunction, damage	①Check if it works properly after resetting ②Check if the memory card is working properly after inserting and removing it ③Need maintenance
Memory card cannot be read	Bad contact or damaged memory card	①Replug the memory card ②Replace the memory card with a new one

Note: It is normal for the App to exit after plugging

06/ Trademark And Legal Statements

“Beaverlab” is a trademark applied for or registered in mainland China by Beaverlabtech Limited Liability Company. and is used on this merchandise. Without the permission of the trademark owner, no person or organization may use the above trademark on goods that are not approved. This manual is produced and copyrighted by Beaverlabtech Limited Liability Company and may not be reproduced or distributed in whole or in any part by any organization or individual without permission. Due to the continuous improvement of product functions, design changes and other reasons, this manual may not match the product you purchased. Please refer to the actual product.

FCC WARNING

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

To maintain compliance with FCC's RF exposure guidelines, This equipment should be installed and operated with minimum 20cm distance between the radiator and your body. Use only the supplied antenna.

Battery Precautions

1. The product uses rechargeable batteries.
2. Rechargeable batteries must be charged under adult supervision.
3. Please use the recommended charger.
4. Power supply electronics shall not short circuit

Microscope Global Parameters

Product Brand	Beaverlab
Product Name	Digital Microscope
Product Model	DDL-M2C
Product Color	White/Gray
Single Cell Capacity	2500mAh Li-ion battery
Material	ABS, PC, Aluminum alloy, LCD, Electronic component

By scanning you will be able to
Download APP
Get exclusive benefits
Unlock more ways to play

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Brand Company: Beaverlabtech Limited Liability Company